


THUASNE

Rebel® Lite
Rebel®
Rebel® Pro

fr	Genouillère ligamentaire rigide articulée	5
en	Hinged rigid ligament knee brace	5
de	Starre Knieführungsorthese	6
nl	Rigide scharnierende ligamentaire kniebrace	7
it	Ginocchiera per legamenti rigida articolata	7
es	Rodillera ligamentaria rígida articulada	8
pt	Joelheira ligamentar rígida articulada	9
da	Stift artikuleret knæbind til ledbånd	9
fi	Nivelletty polven jäykkä nivelsidetuki	10
sv	Ledat styvt knäskydd för ledband	11
el	Επιγονατίδα συνδέσμων άκαμπτη αρθρωτή	11
cs	Pevná kolenní vazová ortéza s kloubem	12
pl	Sztywna, przegubowa orteza stawu kolanowego	13
lv	Stingra ceļa locītavas saišu ortoze ar eņģēm	14
lt	Standus raiščių antkelis su jungtimis	14
et	Liigenduvate küljetugevdustega põlvetugi	15
sl	Pregibna toga opornica za koleno	15
sk	Pevná pohyblivá ortéza na väzy kolenného kĺbu	16
hu	Ízületes merev térdortézis a térdszalagok rögzítésére	17
bg	Съчленена накованка за разтегнати връзки на коляното	17
ro	Genunchieră ligamentară rigidă cu articulație	18
ru	Жесткий лигаментарный ортез для коленного сустава с боковыми шарнирами	19
hr	Čvrsta zglobna ortoza za ligamente koljena	20
zh	硬质关节韧带护膝	20
ar	داعم مفصلي صلب للركبة وأربطتها	21

			
	(a)	(b)	(c)
XS	31 - 39 cm	7,5 - 9 cm	28 - 32 cm
SM	39 - 47 cm	9 - 10 cm	31 - 35 cm
MD	47 - 53 cm	10 - 11,5 cm	34 - 38 cm
LG	53 - 60 cm	11,5 - 12,5 cm	37 - 40 cm
XL	60 - 64 cm	12,5 - 14 cm	38 - 43 cm
XXL	64 - 71 cm	14 - 15,5 cm	43 - 48 cm
XXXL*	71 - 79 cm	15,5 - 16,5 cm	48 - 53 cm

fr	*XXXL disponible uniquement pour la Rebel® Pro.
en	*XXXL available on Rebel® Pro only.
de	*XXXL nur für Rebel® Pro erhältlich.
nl	*XXXL enkel beschikbaar voor de Rebel® Pro.
it	*XXXL disponibile solo per Rebel® Pro.
es	*XXXL disponible solo para Rebel® Pro.
pt	*XXXL disponível somente para a Rebel® Pro.
da	*XXXL findes kun i versionen Rebel® Pro.
fi	*XXXL saatavissa ainoastaan malliin Rebel® Pro.
sv	*XXXL finns endast att tillgå för Rebel® Pro.
el	*XXXL διαθέσιμο μόνο για το μοντέλο Rebel® Pro.
cs	*XXXL k dispozici pouze u Rebel® Pro.
pl	*XXXL dostępne tylko w wersji Rebel® Pro.
lv	*XXXL pieejams tikai Rebel® Pro.
lt	*XXXL dydis galimas tik Rebel® Pro.
et	*XXXL suuruses saadaval ainult Rebel® Pro.
sl	*XXXL na voljo samo Rebel® Pro.
sk	*XXXL k dispozícii len v prípade Rebel® Pro.
hu	Az *XXXL méret csak a Rebel® Pro esetében áll rendelkezésre.
bg	*XXXL се предлага само за Rebel® Pro.
ro	*XXXL disponibilă numai pentru Rebel® Pro.
ru	*Размер XXXL доступен только для модели Rebel® Pro.
hr	*Veličina XXXL dostupna je samo za Rebel® Pro.
zh	*XXXL号仅限Rebel® Pro.
ar	.Rebel® Pro متوفر فقط للمنتج XXXL*

fr	Stabilisation	Contrôle du mouvement
en	Stabilisation	Motion control
de	Stabilisierung	Bewegungsführung
nl	Stabilisatie	Bewegingscontrole
it	Stabilizzazione	Controllo del movimento
es	Estabilización	Control del movimiento
pt	Estabilização	Controlo do movimento
da	Stabilisering	Kontrol af bevægelsen
fi	Stabilioi	Liikkeen hallinta
sv	Stabilisering	Rörelsekontroll
el	Σταθεροποίηση	Έλεγχος της κίνησης
cs	Stabilizace	Kontrola pohybu
pl	Stabilizacja	Kontrola ruchu
lv	Stabilizēšana	Kustības kontrole
lt	Stabilizavimas	Judesio kontrolė
et	Stabiliseerimine	Liikuvuse kontroll
sl	Stabiliziranje	Nadzor gibanja
sk	Stabilizácia	Vymedzenie rozsahu pohybu
hu	Stabilizálás	Mozgásszabályozás
bg	Стабилизация	Контрол на движението
ro	Stabilizare	Controlul mișcării
ru	Стабилизация	Контроль движения
hr	Stabilizacija	Kontrola pokreta
zh	稳定	运动控制
ar	الثابت	التحكم في الحركة

fr	(a) Circonférence de la cuisse (15 cm au-dessus du genou)	(b) Largeur du genou	(c) Circonférence du mollet (15 cm au-dessous du genou)
en	(a) Thigh circumference (15 cm above the knee)	(b) Knee width	(c) Calf circumference (15 cm below the knee)
de	(a) Oberschenkelumfang (15 cm über dem Knie)	(b) Kniebreite	(c) Wadenumfang (15 cm unter dem Knie)
nl	(a) Omvang van de dij (15 cm boven de knie)	(b) Breedte van de knie	(c) Omvang van de kuit (15 cm boven de knie)
it	(a) Circonferenza della coscia (15 cm al di sopra del ginocchio)	(b) Larghezza del ginocchio	(c) Circonferenza del polpaccio (15 cm al di sotto del ginocchio)
es	(a) Circunferencia del muslo (15 cm por encima de la rodilla)	(b) Ancho de la rodilla	(c) Circunferencia de la pantorrilla (15 cm por debajo de la rodilla)
pt	(a) Circunferência da coxa (15 cm acima do joelho)	(b) Largura do joelho	(c) Circunferência da barriga da perna (15 cm abaixo do joelho)
da	(a) Lårets omkreds (15 cm over knæet)	(b) Knæets bredde	(c) Læggens omkreds (15 cm under knæet)
fi	(a) Reiden ympäryys (15 cm polven yläpuolelta)	(b) Polven leveys	(c) Pohkeen ympäryys (15 cm polven alapuolelta)
sv	(a) Omkrets runt låret (15 cm ovanför knäet)	(b) Knäets bredd	(c) Omkrets runt vad (15 cm under knäet)
el	(a) Περίμετρος του μηρού (15 εκ. πάνω από το γόνατο)	(b) Περίμετρος του γόνατος	(c) Περίμετρος της κνήμης (15 εκ. κάτω από το γόνατο)
cs	(a) Obvod stehna (15 cm nad kolénem)	(b) Obvod kolena	(c) Obvod lýtky (15 cm pod kolénem)
pl	(a) Obwód uda (15 cm powyżej kolana)	(b) Szerokość kolana	(c) Obwód łydki (15 cm poniżej kolana)
lv	(a) Augšstilba apkārtmērs (15 cm virs ceļgala)	(b) Ceļgala platums	(c) Apakšstilba apkārtmērs (15 cm zem ceļgala)
lt	(a) Šlaunies apimtis (15 cm virš kelio)	(b) Kelio plotis	(c) Blauzdos apimtis (15 cm žemiau kelio)
et	(a) Reie ümbermõõt (15 cm põlvest kõrgemal)	(b) Põlve ümbermõõt	(c) Sääremarja ümbermõõt (15 cm põlvest allpool)
sl	(a) Obseg stegna (15 cm nad kolonom)	(b) Širina kolena	(c) Obseg goleni (15 cm pod kolonom)
sk	(a) Obvod stehna (15 cm nad kolenom)	(b) Šírka kolena	(c) Obvod lýtky (15 cm pod kolenom)
hu	(a) Comb körfogata (a térd felett 15 cm-rel)	(b) Térd szélessége	(c) Vádli körfogata (a térd alatt 15 cm-rel)
bg	(a) Обиколка на бедрото (15 cm над коляното)	(b) Ширина на коляното	(c) Обиколка на прасеца (15 cm под коляното)
ro	(a) Circumferința coapsei (la 15 cm deasupra genunchiului)	(b) Lățimea genunchiului	(c) Circumferința gambei (la 15 cm sub genunchi)
ru	(a) диаметр окружности бедра (на 15 см выше колена)	(b) Ширина колена	(c) Обхват голени (на 15 см ниже колена)
hr	(a) Opseg bedra (15 cm iznad koljena)	(b) Širina koljena	(c) Opseg lista (15 cm ispod koljena)
zh	(a) - 大腿周长 (膝盖上方15厘米处)	(b)膝关节宽度	(c) 小腿周长 (膝盖下方15厘米处)
ar	(a) محيط الفخذ (15 سم فوق الركبة)	(b) عرض الركبة	(c) محيط ربة الساق (15 سم فوق الركبة)

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HINGED RIGID LIGAMENT KNEE BRACE

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Product available in left or right version.

Available in 3 versions: Rebel® Lite, Rebel®, Rebel® Pro.

Available in 6 sizes for Rebel® Lite and Rebel® and 7 sizes for Rebel® Pro. Rebel® Lite is available in length 33cm, Rebel® in length 33 cm and 36 cm and Rebel® Pro in length 36 cm (other lengths are available in option).

The device is composed of:

- a frame ensuring the brace's rigidity,
- the TM5+ hinge ① mimicking the knee's natural movement, protected by a cover,
- 2 condylar pads,
- 4 posterior half-straps ②, ③, ④, ⑤ (including a synergistic suspension strap ②), 1 anterior half-strap ⑥ and 1 optional anterior half-strap ⑦.

A kit with 2 additional condylar pads, extension limitations ⑧ and a screwdriver is included in the brace's box.

The kit containing flexion limitations ⑨ could be present in the box or can be ordered as an option depending of the version.

Composition

Textile components: elastane - polyamide - ethylene vinyl acetate - silicone.

Rigid components: aluminium - stainless steel - polyoxymethylene - vinyl - polypropylene - polyamide - epoxy resin - brass.

Properties/Mode of action

Stabilization of knee joint ligaments provided by the rigid hinged side uprights.

- Extension can be adjusted to 0°, 5°, 10°, 15°, 20° and 30°.
- Flexion can be adjusted to 0°, 30°, 45°, 60°, 75°, 90° and 110°.

The brace is maintained on the leg thanks to:

- the Synergistic Suspension Strap ②;
- the shape of the tibial bolster (rounded lateral part, flat medial part) that prevents rotation,
- the Compression/Suspension pack providing improved hold on the thigh (optional).

Indications

Conservative treatment of knee ligament injuries and/or ruptures (cruciate and/or lateral ligaments).

Post-operative immobilisation or/and rehabilitation.

Post-traumatic immobilisation.

Joint instability/laxity (including for knee osteoarthritis).

Contraindications

Do not use the product if the diagnosis has not been confirmed.

Do not use in the event of known allergy to any of the components.

Do not apply the product in direct contact with broken skin or an open wound without an adequate dressing.

History of venous or lymphatic disorders.

Do not use in the event of major venous thromboembolic history without thromboprophylaxis.

Precautions

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

In the event of any modification in the product's performance, remove it and consult a healthcare professional.

Before any sports activity, check the compatibility of the use of this medical device with your healthcare professional.

Do not wear the product in a medical imaging machine.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

Undesirable side-effects

This device can cause skin reactions (redness, itching, burns, blisters,etc.) or wounds of various degrees of severity.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

It is recommended to wear the device directly on the skin, unless contraindicated.

Contouring the brace:

To make the knee brace more comfortable and improve its hold, the uprights can be shaped. Apply a slight pressure in the desired direction. Repeat the process on the other upright if needed.

If subsequent shaping of the uprights is needed, consult a healthcare professional.

The use of thicker condylar pads, provided with the brace, may enhance the fit and provide more proprioceptive pressure to the sides of the knee.

Adjusting the strap length:

If one or more of the knee brace straps is too long, each strap can be cut to the desired length.

To do this, remove the self-fastening hook tab from the end of the strap, cut the strap to the desired length, and reattach the self-fastening hook tab to the end of the strap **Ⓢ**.

Take care to not cut any strap too short.

Additionally, any comfort pad that is attached to the inside of a strap must be removed and cut shorter if the pad interferes with your ability to completely tighten the strap.

Flexion/extension control:

The adjustment of the flexion/extension must be decided and performed by the healthcare professional, not by the patient.

Adjusting the extension limitation Ⓢ:

The brace's default position is limitation of extension to 0°.

To change this setting, please follow the following instructions. **to be repeated with the same limitation on both hinges.**

Extension can be limited to 0°, 5°, 10°, 15°, 20° and 30°.

- Choose the desired extension limitation on the corresponding plastic piece.
- Remove the screw located on the side of each hinge.
- Flex the hinge slightly and remove the existing extension limitation. Take care to note the direction this limitation is facing.
- Insert the desired limitations, hole end first, with the hook end at the top and facing forward.

Straightens the hinge to full extension to confirm the extension limitations are correctly positioned. The small hole in each limitation must be aligned and visible through the screw hole so the screw will thread into the limitation.

5. Reinsert and tighten the screw.

Do a few flexions/extensions to make sure that the limitation is properly locked at the desired angle.

Adjusting the flexion limitation Ⓢ:

Flexion can be limited to 0°, 30°, 45°, 60°, 75°, 90° and 110°.

- Choose the desired flexion limitation on the corresponding metallic piece. Each limitation has the degree etched into the surface.
- Remove the two screws from the posterior aspect of each hinge cover and take out the spacer that was installed at the factory. Insert the limitation with the flat end facing up and position it so that both holes in the limitation are visible through the screw holes in the cap. Thread and tighten the screws through the cap and into both holes in the limitation.

Do a few flexions/extensions to make sure that the limitation is properly locked at the desired angle.

Positioning the brace on the leg:

-Have the patient sit on the edge of a chair, and have the patient bend his/ her knee to a position of approximately 30-60° of flexion.

Ensure that the four straps that cross the back of the leg are open **ⓐ**, **ⓑ**, **ⓒ**, **ⓓ**, and to ease application, fold the self-fastening straps back.

-Place the knee brace on the leg.

Ensure that the brace size is appropriate by checking that there is no space between the brace and the leg. If there is any gapping, or if the brace is too tight, if possible, conform the frame of the brace or retake measurements of the leg and select the appropriate size (see size table).

-Identifying knee center:
The condylar pads attached to the inside of the hinges should press against the sides of your knee.

The center of the hinges against the sides of the knee must be aligned with the upper third of the patella.
NOTE: It is always better to put your brace on a little too high than too low.
-First tighten and fasten the synergistic suspension strap **ⓐ** located above the calf.

It is the most important strap for preventing the brace from slipping down your leg. By tightening the synergistic suspension strap **ⓐ** above the calf muscle in the flexion fold behind the knee you will ensure that the hinges are at the correct level on the leg.

-Next fasten the anterior tibial strap **ⓑ** but not too tight.

The anterior-posterior position of the hinges can be modulated by adjusting straps **ⓐ** and **ⓑ**.

The hinges and uprights should be aligned with the midline of the sides of the leg, or just posterior to midline (2nd third of the leg in the sagittal plane).

After you pre-set the anterior strap, instruct the patient to leave it closed when putting on or taking off the brace.
-Next fasten the lower strap **ⓒ**, and then the two straps behind the thigh **ⓓ** and **ⓔ**, and finish by fastening the anterior thigh strap **ⓕ** (optional).
The chafes that you feel the strap through can be angled to ensure the straps fit flat against the back side of your leg.

-As a final step, push down on the brace to make sure it is strapped securely and stays in place.
We also recommend to stretch out the leg and retighten the lower thigh strap **ⓐ** to achieve maximum control.

If the patient feels the need to tighten the straps during an activity, make sure to tighten them following the steps indicated in the fitting instructions.

Care

Close the self-fastening tabs before washing. Hand wash. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

Hinges Ⓢ:

The hinges on the brace are pre-lubricated in the factory. If sand, dirt or water gets inside the hinges, they may require lubricating again. If you notice the hinges not gliding smoothly, a few drops of a synthetic lubricant can be applied. Wipe off any excess lubricant before wearing the brace to prevent stains on clothing.

Straps:

After prolonged use, if the fibers on your strap do not adhere as well to the self-fastening tab, cut the strap shorter so the self-fastening tab adheres to a section of the strap that has fresher fibers. If this is not possible, you should contact the medical provider who fit your brace.

Pads:

The brace is lined with padding that provides a comfortable interface between the leg and the shells. The straps also have pads. Do not remove these pads from the brace or straps. Wipe the pads after each use to remove any moisture and let the pads air dry. You can also clean the pads with a mild anti-bacterial soap and rinse them off with fresh water. Do not wash the pads in a machine and do not dry in a tumble dryer.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

COMMERCIAL WARRANTY AND WARRANTY LIMITATIONS

Thuasne offers a free, limited commercial warranty to the user, in the territory where the device was purchased, against defects in manufacturing and workmanship for a period of:

-six months for the textile components;

-one year for the frame and the brace hinge.

The limited warranty is effective from the date of purchase of the product by the end-user.

The limited commercial warranty does not apply to any defects in manufacturing and workmanship in case of:

-misuse of the product or any damage occurred by a usage outside the normal and intended use of the product as mentioned in the instructions for use,

-damages occurred as part of attempts to modify the product.

Any deterioration or improper trimming of the product during its modification or adjustment process by the healthcare professional when fitting device is expressly excluded from this warranty.

Any claim for this commercial warranty must be sent by the user to the entity where the product was purchased, which will forward this claim to the corresponding Thuasne entity.

Any warranty claim will first be reviewed by Thuasne to determine if the conditions of the limited warranty are fulfilled and do not fall into one of the cases of exclusion of the commercial warranty.

To benefit from the warranty, the buyer must mandatorily provide an original and dated proof of purchase of the product.

If the conditions of the limited warranty are fulfilled and the claim is made by the user or its legal representative (parents, guardian...) within the warranty delays indicated above, the buyer will get a new substitution product.

It is expressly agreed that this commercial warranty is in addition to the legal warranties binding the entity which sold the product to the user, in accordance with the applicable local legislation in the country of purchase of the product.

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©Thuasne - 2031001 (2022-02)