



LombaFirst


Standard lumbar belt

Elastic fabric

The **beaded selvedge of the fabric** ensure a good **comfort** and a certain elasticity at **the ribs and hips areas**

Easy to use

Finger-loops to help easily tighten the belt



Sense of wear


To indicate the **sense of wear** and **facilitate the fitting**, the product have a **logo** on the abdominal panel

Discreet & lightweight

Thin fabric and **seamless** sewing technology

Black color

A **sober color** to wear the brace **discreetly**



Integrated spring-loaded steel back stays

The stays **can be conformed** to the patient lordosis

Smart angles

Smart and ergonomic Design which enable the proper comfort in seating position

Machine wash at 30°

Medical indications

Common low back pain (acute, subacute and chronic phrases), Spinal stenosis

Height

Available in 21 & 26 cm

Sizing

Black



T1	60 - 80 cm
T2	75 - 95 cm
T3	90 - 110 cm
T4	105 - 125 cm
T5	120 - 141 cm

DESIGNATION	EAN13
LOMBAFIRST H21 T1	3111790876106
LOMBAFIRST H21 T2	3111790876113
LOMBAFIRST H21 T3	3111790876120
LOMBAFIRST H21 T4	3111790876137
LOMBAFIRST H21 T5	3111790876144
LOMBAFIRST H26 T1	3111790875109
LOMBAFIRST H26 T2	3111790875116
LOMBAFIRST H26 T3	3111790875123
LOMBAFIRST H26 T4	3111790875130
LOMBAFIRST H26 T5	3111790875147