



LombaStart

Standard Lumbar belt

Elastic fabric

The **beaded selvedge of the fabric ensures** a good **comfort** and a certain elasticity at **the ribs and hips areas**

Easy to use

Finger-loops to help easily tighten the belt



Sense of wear

To indicate the **sense of wear** and **facilitate the fitting**, the product have a **logo** on the abdominal panel

Additional strap

An **immediate relief** of the lumbar area can be felt while tighten the belt

Discreet & lightweight


Thin fabric and seamless sewing technology

Integrated spring-loaded steel back stays

The stays **can be conformed** to the patient lordosis

Black color

A **sober color** to wear the brace **discreetly**

 Machine wash at 30°

Smart angles

Smart and ergonomic Design which enable the proper comfort in seating position



Medical indications

Common low back pain (acute, subacute and chronic phrases), Spinal stenosis

Height

Available in 21 & 26 cm

Sizing

Black



T1	60 - 80 cm
T2	75 - 95 cm
T3	90 - 110 cm
T4	105 - 125 cm
T5	120 - 141 cm

DESIGNATION	EAN13
LOMBASTART 26 T1 NOIR	3111790873112
LOMBASTART 26 T2 NOIR	3111790873129
LOMBASTART 26 T3 NOIR	3111790873136
LOMBASTART 26 T4 NOIR	3111790873143
LOMBASTART 26 T5 NOIR	3111790873150
LOMBASTART 21 T1 NOIR	3111790874119
LOMBASTART 21 T2 NOIR	3111790874126
LOMBASTART 21 T3 NOIR	3111790874133
LOMBASTART 21 T4 NOIR	3111790874140
LOMBASTART 21 T5 NOIR	3111790874157