



**THUASNE**

# Rebel PRO Ligament

Patients who have experienced ligament injuries and/or reconstructive surgery. For rehabilitative or prophylactic use.

## Features


- Optimised & customised fit - Custom assembled to 3 leg measurements
- Enhanced brace positioning via Townsend Motion TM5+Hinge- tracks with the anatomical roll-and-glide motion of the knee to keep the brace properly positioned on the leg.
- Superior suspension via synergistic suspension strap
- Improved suspension via C/S package
- Rotational control via tibial shell contour
- Suspension augmentation via anti-migration straps
- Easy donning via clip buckles
- Suitable for heavy users and contact sports- rigid aircraft aluminium frame

## Options


- Extension assist kit- suitable for patients who present with quadriceps weakness or fatiguable quadriceps
- Frame height options
- Frame configuration options
- Flexion Stop Kit
- Range of colours



## Anti-migration Technology



**C/S Package**  
(Compression/Suspension)



**Synergistic Suspension Strap**

## Powdercoat Colours



To order a Rebel Pro please simply complete form (overleaf)



# Rebel Pro ORDER FORM

### Account Contact Information

Clinician Name: \_\_\_\_\_ Clinician Email: \_\_\_\_\_ Clinician Phone: \_\_\_\_\_  
Account Name: \_\_\_\_\_ Purchase Order Number: \_\_\_\_\_

### For Thuasne Office Use Only

Account Number: \_\_\_\_\_ Order Number: \_\_\_\_\_  
Purchase Order Number: \_\_\_\_\_ Purchase Order Number: \_\_\_\_\_

### Patient Information

Fit Date: \_\_\_\_\_

Patient's Last Name: \_\_\_\_\_

Patient's First Name: \_\_\_\_\_

Male  Female Age \_\_\_\_\_

Weight \_\_\_\_\_ (LBS) Height \_\_\_\_\_ (IN)

Leg:  Left  Right

### Brace Configuration

Thigh Shell Length  7" anterior  8" anterior

Tibia Shell Length  6" anterior  7" anterior  
 8" anterior  7" posterior

### Colour

#### Matte Finish

Black (Standard)  Grey  White  
 Atlantic (Light Blue)  Red

#### Satin Finish

Lemon  Orange  Lime  
 Fuchsia  Pacific (Dark Blue)

### Accessories

- Flexion Stop Kit (Extension stops included) £10 charge
- Extension assist bands/posts\* £75 charge
- Quick release buckles
- Anti-migration silicon infused strap pads
- Spooner patella stabilizing attachment\* £65 charge
- C/S Package

### Sizing

If your patient has proportional leg sizing (see sizing reference, below) or if you are ordering a brace for stock inventory please select from the size options.

X-Small  Small  Medium  Large  X-Large

Size	6" Above	M-L Knee Width	6" Below
XS	12.5 - 15.5"	3 - 3.5"	11 - 12.5"
SM	15.5 - 18.5"	3.5 - 4"	12.25 - 13.25"
MD	18.5 - 21"	4 - 4.5"	13.25 - 15"
LG	21 - 23.5"	4.5 - 5"	14.25 - 15.75"
XL	23.5 - 25"	5 - 5.5"	15 - 17"
XXL	25 - 28"	5.5 - 6"	17 - 19"
XXXL*	28 - 31"	6 - 6.5"	19 - 21"

\*XXXL available on Rebel Pro only

### Measurement Data

#### Required for "custom-fitted" Assembly\*

Custom-fitted (Fabricated from three leg measurements)\*


If your patient has a disproportionate leg, Thuasne USA will customize the assembly of your patient's brace. Please provide leg measurements beside the illustration, below.

More than one size difference between thigh and calf may require a full custom knee brace - please refer to the SpryStep Knee Brace

\_\_\_\_\_ Circumference 6 inches above mid-patella

\_\_\_\_\_ Medial-Lateral Knee Width (not circumference) at knee center

\_\_\_\_\_ Circumference 6 inches below mid-patella



Special Instructions: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

Thuasne UK Ltd

Unit 4 Orchard Business Centre, North Farm Road, Tunbridge Wells, TN2 3XF

Tel: 01892 481620 - Fax: 01892 481609 - Email: orders@thuasne.co.uk

www.thuasne.co.uk