



THUASNE

Rebel Standard Ligament


Patients who have experienced ligament injuries and/or reconstructive surgery. For rehabilitative or prophylactic use.

Features


- Stock sizes see chart
- Enhanced brace positioning via Townsend Motion TM5+ Hinge- tracks with the anatomical roll-and-glide motion of the knee to keep the brace properly positioned on the leg.
- Superior suspension via synergistic suspension strap
- Improved suspension via C/S package
- Rotational control via tibial shell contour
- Suspension augmentation via anti-migration straps
- Easy donning via clip buckles
- Suitable for the majority of patients
- ACL configuration



Anti-migration Technology



C/S Package
(Compression/Suspension)



Synergistic Suspension Strap

Size	6" Above	M-L Knee Width	6" Below
XS	12.5 - 15.5"	3 - 3.5"	11 - 12.5"
SM	15.5 - 18.5"	3.5 - 4"	12.25 - 13.25"
MD	18.5 - 21"	4 - 4.5"	13.25 - 15"
LG	21 - 23.5"	4.5 - 5"	14.25 - 15.75"
XL	23.5 - 25"	5 - 5.5"	15 - 17"
XXL	25 - 28"	5.5 - 6"	17 - 19"
XXXL*	28 - 31"	6 - 6.5"	19 - 21"

*XXXL available on Rebel Pro only. See Rebel Pro datasheet

Thuasne UK Ltd

Unit 4 Orchard Business Centre, North Farm Road, Tunbridge Wells, TN2 3XF
Tel: 01892 481620 - Fax: 01892 481609 - Email: orders@thuasne.co.uk
www.thuasne.co.uk