



THUASNE

New Design

GenuPro Comfort

**ELASTIC PATELLAR
KNEE SUPPORT**

Comfortable⁽¹⁾
patella alignment
and joint support



Anterior knee pain, patellar discomfort,
mild osteoarthritis

GenuPro Comfort

ELASTIC PATELLAR KNEE SUPPORT

- Soft, comfortable⁽¹⁾ and breathable⁽²⁾
- Maintains and guides patella⁽³⁾
- Provides joint support⁽³⁾


Soft & comfortable

Pillow-soft viscose yarn

Also called “synthetic silk”, the qualities of the viscose yarn provide the brace’s softness and comfort⁽⁴⁾

Specific knit at the back of the knee

Stretchable knit in popliteal area enhances comfort especially during flexion movements and facilitates application and removal





Patella support

AirFlex insert

- Maintains patella's alignment thanks to oval-shaped design⁽³⁾
- Conforms to the patella via flexible structure with triangular-shaped pattern
- Allows air to flow through the openings and grooves
- Thicker at the top



Compression

3D compressive knit

- Reinforces awareness and secures movement through proprioceptive effect
- Helps contain and reduce edema⁽³⁾
- Stabilizes⁽³⁾ the joint in combination with lateral flexible stays



3D Pull-tabs

Facilitates grip and application of brace



PATIENT PROFILE

- Patients requiring use on a regular basis
- Patients who value comfort while wearing


New Design

GenuPro Comfort

Soft, comfortable⁽¹⁾ and breathable⁽²⁾


Indications

- Anterior knee or patella pain
- Mild patellar instability
- Patella tendinopathy (including jumper's knee syndrome)
- Knee pain and/or swelling management
- Knee osteoarthritis (mild to moderate)



Sizes

Available in 7 sizes, reflecting current morphological trends



	(a)	(b)	ref
0	35 - 37 cm	29,5 - 32 cm	234604 20009899
1	37 - 41 cm	32 - 34,5 cm	234604 20109899
2	41 - 45 cm	34,5 - 37 cm	234604 20209899
3	45 - 50 cm	37 - 39,5 cm	234604 20309899
4	50 - 55 cm	39,5 - 42 cm	234604 20409899
5	55 - 61 cm	42 - 44,5 cm	234604 20509899
6	61 - 67 cm	44,5 - 47 cm	234604 20609899



www.thuasne.com

(1) Thanks to viscose yarn, anatomical shaped knit, specific knitted area at the back of the brace and popliteal area
(2) IN2103 - Evaluation RET (July 2021)
(3) Internal CE marking data
(4) Soft & Care : Techni sens EH2798 study for Thuasne - February 2019. Ease of application and removal, comfort, softness and purchase intention for the Thuasne Soft&Care class 3 product

GenuPro Comfort availability might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region. Read the product's instructions for use, indications and contraindications carefully. The Medical Device mentioned on this document is CE marked according to the European council regulation 2017/745. ref. : 2201006 Last revision date: 01/2022

Photos : Studio Caterin, Taves/shutterstock.com

