



THUASNE

## SpryStep® OA Knee

### Contact Information

☐ Clinician ☐ Fitter/Assistant/Tech ☐ Other: \_\_\_\_\_

Name: \_\_\_\_\_

Email: \_\_\_\_\_ Phone: \_\_\_\_\_

### Ordering Clinician

☐ CPO ☐ CO ☐ CP ☐ Other: \_\_\_\_\_

Name: \_\_\_\_\_

Email: \_\_\_\_\_ Phone: \_\_\_\_\_

### Billing & Shipping

PO#: \_\_\_\_\_

Billing Account#: \_\_\_\_\_

Shipping Account#: \_\_\_\_\_

Shipping Address: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

### Shipping Preference

☐ Ground

☐ Next Day AM

☐ Next Day PM

☐ 2-Day AM

☐ 2-Day PM

(If no preference is indicated, this order will be shipped 2 Day PM) Note: We do not ship products directly to patients.

### Patient Information

Fit Date: \_\_\_\_\_ Patient ID: \_\_\_\_\_

Age \_\_\_\_\_ ☐ Male ☐ Female

Weight \_\_\_\_\_ ☐ Lbs. ☐ Kg. Height \_\_\_\_\_ ☐ in. ☐ cm.

Leg: ☐ Left ☐ Right

Diagnosis: (ex: Ligament laxity, ROM limitations, etc.) \_\_\_\_\_

Compartment: ☐ Unload Medial ☐ Unload Lateral


### Measurement Data

These measurements are required to check the accuracy of the patient model submitted, a patient model must be provided for fabrication (scan).

Proximal circumference  
7 in / 175mm above mid-patella

Medial-Lateral Knee Width  
(not circumference) at knee center

Distal circumference  
7 in / 175mm below mid-patella



### Brace Configuration

NB by default: riveted anchor tabs + d-rings + 1/4" padding + condylar pads + 2 additional thicker condylar pads + synergistic suspension strap

Hinge (Extension stop kit included with hinges)

- ☐ TM5 Aluminum
- ☐ TM5 Aluminum without Loadshifter, Correction built into the frame: \_\_\_\_\_° ☐ Varus ☐ Valgus
- ☐ TM5 Stainless without Loadshifter, Correction built into the frame: \_\_\_\_\_° ☐ Varus ☐ Valgus

### Optional Hinge Accessories

- ☐ Flexion stop kit\* ☐ Extension assist bands/posts\*

### Brace Rigidity/Stiffness

For larger and heavier framed patients- increased rigidity/stiffness is recommended

- ☐ Level 1 (default) ☐ Level 2 (medium) ☐ Level 3 (high)  
recommended for Impact Sports

### Femoral shell length

- ☐ 7 in (default) 175mm ☐ -1 in 150mm ☐ +1 in 200mm ☐ +2 in 225mm

☐ Other \_\_\_\_\_

### Femoral shell configuration

☐ Anterior



☐ Posterior



### Tibial shell length

- ☐ 7 in (default) 175mm ☐ -1 in 150mm ☐ +1 in 200mm ☐ +2 in 225mm

☐ Other \_\_\_\_\_

### Tibial shell configuration

☐ Anterior



☐ Posterior



### Options

- ☐ Full Shell\*
- ☐ C/S Package (not available with FullShell)
- ☐ Combined Instability Strap (PCL)

### Accessories

- ☐ Spooner Patella Stabilizing Attachment\*
- ☐ Brace cover (pull-on style)\*
- ☐ Cotton Undersleeve 18 in (46cm)\*
- ☐ Neoprene Undersleeve 18 in (46cm)\*
- ☐ Neoprene Undersleeve 22 in (56cm)\*
- ☐ CS wrap\*



Full Shell

\*Indicates additional charges apply

Received Date: \_\_\_\_\_ Thuasne USA's shipping department use only

By filling this order form and placing an order for this device, I hereby certify that I am authorized to dispense this medical device in virtue of any national law governing the fitting and adjustment of orthopedic medical devices.

Please do not provide any personal information (name etc) regarding the patient, but only provide health information necessary to the fabrication of this medical device.

Distributed by Thuasne USA  
4615 Shepard Street, Bakersfield, CA 93313  
Tele: 800.432.3466 - Fax: 800.798.2722  
ThuasneUSA.com

OF-000 REV A