



THUASNE

Rebel Lock

Product Code:
Standard - U30202
Semi-Custom - U30201

Indications

- Patients with significant instabilities; weak extensor muscles, failed TKR, limiting or preventing knee flexion.
- A lightweight option for patients who need stability and functionality to control significant instabilities.

Features

- Superior suspension via synergistic suspension strap
- Improved suspension via C/S package
- Rotational control via tibial shell contour
- Suspension augmentation via anti-migration straps
- Lightweight Aircraft Aluminium Frame
- Drop Lock With Free Motion
- Adjustable Locking Position

Options

- Optimised & customised fit - Custom assembled to 3 leg measurements
- Extension assist kit - suitable for patients who present with quadriceps weakness or fatiguable quadriceps
- Frame height options
- Flexion stops at 30, 60, or 90 degrees
- Range of colours (stock frame colour is black)



Extension Assist



This optional extension assistance mechanism uses strong latex bands slipped over two posts – one proximal and one distal – to the knee centre. There is no limitation in flexion.

Powdercoat Colours



*For any brace ordered with optional extras please complete form (overleaf)



Rebel Lock ORDER FORM

Please complete all fields to avoid potential delays in processing your order

Account Contact Information

Clinician Name: _____ Clinician Email: _____ Clinician Phone: _____
Account Name: _____ Purchase Order Number: _____

For Thuasne Office Use Only

Account Number: _____ Order Number: _____
Purchase Order Number: _____ Purchase Order Number: _____

Patient Information

Fit Date: _____
Patient's Last Name: _____
Patient's First Name: _____

Male Female Age _____
Weight _____ (LBS) Height _____ (IN)
Leg: Left Right

Brace Configuration

Thigh Shell Length 7 Inch 8 Inch
Tibia Shell Length 7 Inch 8 Inch
Locking Position *(if not indicated, zero degree lock position will be applied)*
 Zero Degree Locking Position
 Five Degree Locking Position

Colour

Matte Finish

Black *(Standard)* Grey White
 Atlantic *(Light Blue)* Red

Satin Finish

Lemon Orange Lime
 Fuchsia Pacific *(Dark Blue)*

Accessories

Flexion Stop Kit* *(installed upon request)* £10 charge
 30° 60° 90°
 Extension assist bands/posts* £75 charge
 Anti-migration silicon infused strap pads
 C/S Package

Sizing

If your patient has proportional leg sizing (see sizing reference, below) or if you are ordering a brace for stock inventory please select from the size options.

X-Small Small Medium Large X-Large

Sizing Reference			
	6" above	M-L width	6" below
X-Small	12.5" to 15.5"	3" to 3.5"	11" to 12.5"
Small	15.5" to 18.5"	3.5" to 4"	12.25" to 13.75"
Medium	18.5" to 21"	4" to 4.5"	13.25" to 15"
Large	21" to 23.5"	4.5" to 5"	14.25" to 15.75"
X-Large	23.5" to 25"	5" to 5.5"	15" to 17"


Measurement Data

Required for "custom-fitted" Assembly*

Custom-fitted *(Fabricated from three leg measurements)**

If your patient has a disproportionate leg, Thuasne USA will customize the assembly of your patient's brace. Please provide leg measurements beside the illustration, below.

_____ Circumference 6 inches above mid-patella
_____ Medial-Lateral Knee Width (not circumference) at knee center
_____ Circumference 6 inches below mid-patella



Special Instructions: _____

Thuasne UK Ltd

Unit 4 Orchard Business Centre, North Farm Road, Tunbridge Wells, TN2 3XF
Tel: 01892 481620 - Fax: 01892 481609 - Email: orders@thuasne.co.uk
www.thuasne.co.uk