



THUASNE

Rebel Reliever[®]

Off-loading knee brace

U03309

Indications

- Moderate to severe unicompartement femoral-tibial osteoarthritis of the knee.
- Unicompartement knee joint degeneration requiring offloading.
- Post-operative offloading of a single compartment (articular or meniscus cartilage defect or osteochondral defect).
- Knee joint instability after injury or surgery on the cruciate and/or collateral ligaments (ACL, PCL, LCL, MCL).

Target patient

Arthritic patient requiring maximum offloading

Features

Offloading of the damaged compartment: 3-pressure point system with Loadshifter system

Freedom of movement: TM5+ hinge reproducing the knee's natural movements.*

Open zone at the back of the knee

Hold: suspension strap above the calf and specially shaped tibial shell to prevent rotation




Comfort: inner foam lining

Practical: fitting buckle and tool-less adjustment of offloading level (Snap lock system)

2 models: right and left

Options

Also available in a custom format. Please see order form on overleaf."



STANDARD SIZES	XS	NA	3 - 3.5"	11 - 12.5"
	S	15.5 - 18.5"	3.5 - 4"	12.25 - 13.25"
	M	18.5 - 21"	4 - 4.5"	13.25 - 15"
	L	21 - 23.5"	4.5 - 5"	14.25 - 15.75"
	XL	23.5 - 25"	5 - 5.5"	15 - 17"
	XXL	25 - 28"	5.5 - 6"	17 - 19"

Rebel Reliever

ORDER FORM

LIGAMENT BRACING SOLUTIONS

Please complete all fields to avoid potential delays in processing your order.

Account Contact Information

Clinician Name: _____ Clinician Email: _____ Clinician Phone: _____

Account Name: _____ Purchase Order Number: _____

For Thuasne Office Use Only

Account Number: _____ Order Number: _____

Purchase Order Number: _____ Purchase Order Number: _____

Patient Information

Fit Date: _____

Patient's Last Name: _____

Patient's First Name: _____

Male Female Age _____

Weight _____ (LBS) Height _____ (IN)

Leg: Left Right

Compartment

Medial Compartment Lateral Compartment

Universal (select when ordering stock inventory)
Medial QA pads installed, Lateral QA pads included in kit.

Shell Length Thigh shell extends 7" to 9" on all braces

Tibia 6" 7" 8"

TM5+ Hinges Includes Extension Stops

Optional Flexion Stop Kit*

Colour

Matte Finish

Black (Standard) Grey White

Atlantic (Light Blue) Red

Satin Finish

Lemon Orange Lime

Fuchsia Pacific (Dark Blue)

Accessories

C/S Package

Sizing

If your patient has proportional leg sizing (see sizing reference, below) or if you are ordering a brace for stock inventory please select from the size options.

X-Small Small Medium Large X-Large

Size	6" Above	M-L Knee	6" Below
XS	NA	3 - 3.5"	11 - 12.5"
SM	15.5 - 18.5"	3.5 - 4"	12.25 - 13.25"
MD	18.5 - 21"	4 - 4.5"	13.25 - 15"
LG	21 - 23.5"	4.5 - 5"	14.25 - 15.75"
XL	23.5 - 25"	5 - 5.5"	15 - 17"
XXL	25 - 28"	5.5 - 6"	17 - 19"

Measurement Data

Required for "custom-fitted" Assembly*

Custom-fitted (Fabricated from three leg measurements)*


If your patient has a disproportionate leg, Thuasne USA will customize the assembly of your patient's brace. Please provide leg measurements beside the illustration, below.

More than one size difference between thigh and calf may require a full custom knee brace - please refer to the SpryStep Knee Brace

_____ Circumference 6 inches above mid-patella

_____ Medial-Lateral Knee Width (not circumference) at knee center

_____ Circumference 6 inches below mid-patella



Special Instructions: _____

Thuasne UK Ltd

Unit 4 Orchard Business Centre, North Farm Road, Tunbridge Wells, TN2 3XF

Tel: 01892 481620 - Fax: 01892 481609 - Email: orders@thuasne.co.uk

www.thuasne.co.uk