



THUASNE

**HEALTH EXPERTISE
AT THE ATHLETES
SERVICE**

SPORT



Thuasne Group is 175 years old

From its Saint-Etienne base, THUASNE Group has been drawing on its know-how for more than 175 years to conceive, develop and manufacture medical devices enabling all patients and consumers to become active players in their own health. Derived from the shared history of six generations of entrepreneurs, THUASNE offers concrete, simple, adapted and innovative health solutions for improved mobility in more than 85 countries.

Our expertise is aimed at patients in the fields of orthopedics (belts, braces), medical compression (phlebology, lymphology, wound healing), women's health (external breast forms, lingerie) and homecare.

Thuasne today is...



2400
employees globally



A commercial presence in
85 countries



14
industrial sites in Europe and
the United States



22000
professional customers



30 to 40
new product references
every year



One Thuasne brace sold every
2 seconds
in the world



15 million
kilometers of knitted thread



5700
kilometers of fabric produced per year

A family-owned company, Thuasne puts its know-how and high standards at the service of athletes

More than **80 years of expertise**
in the field of sports



THUASNE®: A range of services

If you are a distributor, Thuasne® can offer to provide you with a whole range of marketing tools to help you in the development of your business and make your life easier!

DEDICATED AREA MANAGER

A qualified area manager expert in your commercial zone will help you in the development of your business.

Market understanding, product analysis & registration, portfolio analysis, training... We can help you achieve your objectives with a wide range of expertise and support!

MARKETING LIBRARY (DATA ASSET MANAGEMENT (DAM))


Gain agility and performance!

Thuasne now provides a dematerialised marketing resources platform¹. You will be able to access everything you need to promote Thuasne® products on your brochures, your corporate website, your marketplace or any other support for you or your customers from the <https://mediacenter.thuasne.com/portal/signin> website:

- HD product images.
- HD packaging images.
- HD lifestyle images.
- Thuasne® catalogues and brochures.
- Product Instructions for Use (IFUs) and Regulatory files.
- Product presentations.
- Product promotional videos.
- Product fitting videos².

Nothing could be easier to access Thuasne's resources:

1. Visit the website:
<https://mediacenter.thuasne.com/portal/signin> and fill in your login information (after registration)
2. Carry out your search using the product name or its code (e.g. for Action Reliever: code 234901)
3. Select the files you wish to download
4. Check the rights of use of the marketing resource and choose your asset format
5. Acknowledge having read the Terms and Conditions of Use and agree to respect them
6. Download your marketing materials





ORDER PORTAL EXPORT

Gain speed and autonomy!

In certain countries, it is possible for you to order Thuasne® products directly from our ordering platform (available at this address: www.partner.thuasne.com).

This new portal allows you to:

- Browse through all our products and categories¹
- Create purchase orders
- Create draft orders
- Define your favourite categories and products
- View your order history
- View or download your shipping documents
- View or download your invoices

DESIGN SUPPORT

Always ready to support you!

For any specific need, our Design Department is at your disposal to support you and to promote our partnership through various marketing tools (brochure, leaflets, POS, communication on tradeshows).

Do not hesitate to express your request at any time and send it to your usual contact.



(1) Access and rights to these platforms will be managed by your Area Manager according to your needs.

(2) Marketing resources will be added progressively and according to their type on the Digital Asset Management (DAM).

CONTENTS

SUPPORT BRACES	8
Knee	9
Ankle	12
Hand, wrist	14
Elbow, Shoulder, Back	15
SPORTS COMPRESSION	18
Runtrail UP Activ	20
UP Recovery	21
Team Sport Range	22
CARE & PROTECTION	24
Support Strapping.....	25
Cold & Heat Therapy.....	26
RESELLER OPPORTUNITIES	27





BRACES AND SUPPORT

ADAPTED TO DIFFERENT ATHLETE'S PATHOLOGIES

Thuasne offers 4 categories of braces specially adapted to sports activities:

- Knee
- Ankle
- Hand, wrist
- Shoulder, elbow, back

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

All the medical devices mentioned on this document are CE marked according to the Regulation 2017/745 on medical devices, unless specifically identified as "NOT CE MARKED".

Some of the products featured in this document are also UKCA-marked according to S.I 2002/618 regulation on medical devices.

Devices mentioned in this document are CE class I. Please contact Thuasne should you need any additional information on devices classification.

Please read carefully the instructions for use, indications and contraindications of the product.




UK Responsible Person (UKRP):

Thuasne UK Ltd
Unit 4 Orchard Business Centre
North Farm Road
Tunbridge Wells TN2 3XF,
United Kingdom



Thuasne
120 rue Marius AUFAN 92300 Levallois - Perret
www.thuasne.com



Thuasne Sport patellar bandage

Patella tendinopathy: Jumper's knee and Runner's knee

032902



INDICATIONS

- Joint or ligament inflammation of the kneecap.
- Patellar tendinitis (Runner's knee and Jumper's knee).
- Mild patellar instability.
- Anterior knee or patella pain.

MODE OF ACTION

The Thuasne Sport patellar bandage maintains the patella and helps dampen vibrations by exerting pressure on the patellar tendon of the knee. The pressure applied to the area supports the patella and helps to reduce pain.

COMPOSITION

Rigid components: polyamide, polyoxymethylene, silicone.

Textile components: polyamide, polyurethane, polyethylene, elastane, synthetic rubber (SBR), polyester.



Universal size

1	26 - 50 cm	3111790047278
---	------------	---------------

IT BAND Thuasne Sport

Tightening band for "Iliotibial band syndrome (ITBS)" 032601



INDICATIONS

- Pain or tendinitis of the iliotibial band (so-called "Iliotibial band syndrome" or "ITB syndrome").
- Sport activities resumption.

MODE OF ACTION

Localized compression on the iliotibial band thanks to the removable insert.

Proprioceptive effect thanks to the combination of the anti-slip and the tightening system providing compression and support.

Adjustment aid thanks to the tightening buckle.

Held in place by the hook and loop system.

COMPOSITION

Rigid components: polyamide, polyoxymethylene.

Textile components: elastane, synthetic rubber (polychloroprene), polyester, thermoplastic polyurethane, ethylene vinyl acetate (EVA).



Universal size

1	27 - 51 cm	3111790046974
---	------------	---------------

Thuasne Sport thigh support

Thigh muscle injuries

057203



UK
CA



Support

INDICATIONS

- Muscular thigh injuries.
- Sport activities resumption.

MODE OF ACTION

Fabric allowing an even compression.

Muscle support.

Limitation of the heat loss by the synthetic rubber (SBR)⁽¹⁾.

COMPOSITION

Synthetic rubber (SBR), polyamide, elastane.



S	44 - 48 cm	3111795722019
M	49 - 53 cm	3111795722026
L	54 - 58 cm	3111795722033
XL	59 - 63 cm	3111795722040

Thuasne Sport open knee brace

Pain or unstable knee

057303



UK
CA



Support

INDICATIONS

- Knee pain and/or swelling management.
- Joint instability/laxity.
- Sport activities resumption.

MODE OF ACTION

Limitation of the heat loss by the synthetic rubber⁽¹⁾.

COMPOSITION

Synthetic rubber (chloroprene), polyamide, elastane.



S	31 - 33 cm	3111795732018
M	34 - 36 cm	3111795732025
L	37 - 39 cm	3111795732032
XL	40 - 42 cm	3111795732049
XXL	43 - 45 cm	3111795732056

(1) Internal CE marking data.

(1) Internal CE marking data.

Thuasne Sport reinforced knitted knee support

Mild sprain, unstable knee

035401



←
This product contains latex

UK
CA



Support

INDICATIONS

- Knee pain and/or swelling management.
- Joint instability/laxity.
- Sport activities resumption.

MODE OF ACTION

Elastic knit ensuring proprioception.

COMPOSITION

Textile components: polyamide, elastane, polyester, elastodiene, viscose, natural latex, silicone.

Rigid components: polyurethane, polyethylene, steel.



S	32 - 34 cm	3111790354000
M	35 - 37 cm	3111790354017
L	38 - 41 cm	3111790354024
XL	42 - 45 cm	3111790354031
XXL	46 - 49 cm	3111790354048

Thuasne Sport reinforced knee support

Pain or unstable knee

057003



UK
CA



Support

INDICATIONS

- Knee pain and/or swelling management.
- Joint instability/laxity.
- Sport activities resumption.

MODE OF ACTION

Knee support with 2 flexible side stays.

Limitation of the heat loss by the synthetic rubber (chloroprene)⁽¹⁾.

COMPOSITION

Rigid component: steel.

Textile components: synthetic rubber (chloroprene), polyamide, elastane, polyester, natural rubber.



S	31 - 33 cm	3111795702011
M	34 - 36 cm	3111795702028
L	37 - 39 cm	3111795702035
XL	40 - 42 cm	3111795702042
XXL	43 - 45 cm	3111795702059

(1) Internal CE marking data.

Thuasne Sport knee support

Knee pain and/or swelling management

033403



UK
CA



Support

INDICATIONS

- Knee pain and/or swelling management.
- Sport activities resumption.

MODE OF ACTION

Joint compression thanks to the elastic knit improving proprioception.

COMPOSITION

Polyamide, elastane, cotton.



S	32 - 34 cm	3111790334316
M	35 - 37 cm	3111790334323
L	38 - 41 cm	3111790334330
XL	42 - 45 cm	3111790334347
XXL	46 - 49 cm	3111790334354

Thuasne Sport reinforced patellar knee brace

Mild tendinopathy or patellar instability

035501



UK
CA



Support

INDICATIONS

- Anterior knee or patella pain.
- Mild patellar instability.
- Patella tendinopathy (including jumper's knee syndrome).
- Knee pain and/or swelling management.
- Sport activities resumption.

MODE OF ACTION

Stabilization of the patella.

Uniform knee compression: anatomically-shaped knit.

Fit on the leg: and anti-slip system, flexible lateral stays.

Easy to fit: comfort zone over the patella.

Comfort zone over the popliteal fossa.

COMPOSITION

Textile components: polyamide, polyester, polyurethane, elastane, silicone, ethylene vinyl acetate.

Rigid components: steel, polyethylene



S	32 - 34 cm	3111790355113
M	35 - 37 cm	3111790355120
L	38 - 41 cm	3111790355137
XL	42 - 45 cm	3111790355144
XXL	46 - 49 cm	3111790355151

Thuasne Sport open strapping knee brace

Pain or unstable knee

033103



← This product contains latex

UK
CA

Support

INDICATIONS

- Joint instability/knee laxity.
- Knee pain and/or swelling management.
- Sport activities resumption.

MODE OF ACTION

Stabilization of the knee ligaments provided by the double strapping system.

Anti-slip system to prevent the knee brace slipping.

Opening to make fitting easier.

Limitation of the heat loss by the synthetic rubber (SBR)⁽¹⁾.

COMPOSITION

Textile components: polyamide, synthetic rubber (SBR), polyester, thermoplastic polyurethane, elastane, ethylene vinyl acetate, polyethylene, natural latex.

Rigid components: polyamide, zinc alloy.



S	32 - 36 cm	3111793313318
M	37 - 40 cm	3111793313325
L	41 - 44 cm	3111793313332
XL	45 - 48 cm	3111793313349

Thuasne Sport Novelastic knee strap

Pain or unstable knee

027002

UK
CA

Support

INDICATIONS

- Joint instability/knee laxity.
- Knee pain and/or swelling management.
- Sport activities resumption.

MODE OF ACTION

Stabilization of the knee ligaments provided by the double strapping system.

An anatomically-shaped knitting for even compression and proprioception.

COMPOSITION

Rigid components: polyamide, steel, polyethylene.

Textile components: polyamide, polyester, polyurethane, elastane, elastodiene, silicone.



S	32 - 34 cm	3111792700317
M	35 - 37 cm	3111792700324
L	38 - 41 cm	3111792700331
XL	42 - 45 cm	3111792700348
XXL	46 - 49 cm	3111792700355

Thuasne Sport reinforced ligament knee brace

Knee ligament injuries

033503



Support

INDICATIONS

- Preventive use for the knee ligament injuries (lateral ligaments).
- Joint instability/knee laxity.
- Knee pain and/or swelling management.
- Sport activities resumption.

MODE OF ACTION

Stabilization of the knee ligaments thanks to the rigid hinged side uprights.

Uniform compression: anatomically-shaped knit.

Staggered straps to prevent tourniquet effect.

Held on the leg by the suspension strap (lower strap).

COMPOSITION

Textile components: polyamide, polyester, elastane, polyurethane, thermoplastic polyurethane, synthetic rubber (SBR), natural rubber.

Rigid components: polyamide, polyoxymethylene.



S	43 - 46 cm	3111790046981
M	46 - 49 cm	3111790046998
L	49 - 52 cm	3111790047001
XL	52 - 55 cm	3111790047018



Thuasne Sport ankle support

Mild sprain or unstable ankle

033303



UK
CA



Support

INDICATIONS

- Joint instability/laxity.
- Proprioceptive support for rehabilitation.
- Mild ankle sprain (acute phase).
- Sport activities resumption.

MODE OF ACTION

Muscle restraint at the ankle and comfort zone around the malleoli.
Uniform compression thanks to the elastic knit that improves proprioception.

COMPOSITION

Polyamide, elastane, cotton.



S	19 - 21 cm	3111790333319
M	21 - 23 cm	3111790333326
L	23 - 25 cm	3111790333333
XL	25 - 27 cm	3111790333340
XXL	27 - 30 cm	3111790333357

(1) Internal CE marking data.



Thuasne Sport open ankle support

Mild sprain and unstable ankle

057103



UK
CA



Support

INDICATIONS

- Joint instability/laxity.
- Proprioceptive support for rehabilitation.
- Mild ankle sprain (acute phase).
- Sport activities resumption.

MODE OF ACTION

Limitation of the heat loss by the synthetic rubber (chloroprene)⁽¹⁾.

COMPOSITION

Synthetic rubber (CR), polyamide, elastane.



S	19 - 20 cm	3111795712010
M	21 - 22 cm	3111795712027
L	23 - 24 cm	3111795712034
XL	25 - 27 cm	3111795712041

Thuasne Sport reinforced ankle support

Mild sprain or unstable ankle

035301



UK
CA



Support

INDICATIONS

- Joint instability/laxity.
- Proprioceptive support for rehabilitation.
- Mild ankle sprain (acute phase).
- Sport activities resumption.

MODE OF ACTION

Ankle compression thanks to the elastic knit.
Anatomical cut with a comfort zone around the malleoli.

COMPOSITION

Polyamide, elastodiene, elastane, viscose.



S	19 - 21 cm	3111790353003
M	21 - 23 cm	3111790353010
L	23 - 25 cm	3111790353027
XL	25 - 27 cm	3111790353034
XXL	27 - 30 cm	3111790353041

Thuasne Sport strapping ankle brace

Mild sprain or unstable ankle

033003

UK
CA

Support

INDICATIONS

- Joint instability/laxity.
- Proprioceptive support for rehabilitation.
- Mild ankle sprain (acute phase).
- Sport activities resumption.

MODE OF ACTION

Uniform ankle compression: elastic knit.

Limitation of extreme lateral movements: double elastic straps.

Promotes the evacuation of perspiration: elastic perforated envelope.

COMPOSITION

Polyamide, synthetic rubber (SBR), polyester, thermoplastic polyurethane, elastane.



S	19 - 22 cm	3111793303319
M	22 - 25 cm	3111793303326
L	25 - 28 cm	3111793303333
XL	28 - 30 cm	3111793303340

Thuasne Sport Novelastic strapping ankle brace

Mild sprain or unstable ankle

028002

UK
CA

Support

INDICATIONS

- Joint instability/laxity.
- Proprioceptive support for rehabilitation.
- Mild ankle sprain (acute phase).
- Sport activities resumption.

MODE OF ACTION

Ankle support: bi-directional "8" strapping to preserve bending/extension.

Ankle compression exerted by the elastic knit to improve proprioception. Anatomical cut and comfort zone around the malleoli.

COMPOSITION

Polyamide, elastodiene, elastane, polyurethane, polyester.



S	19 - 21 cm	3111792800314
M	21 - 23 cm	3111792800321
L	23 - 25 cm	3111792800338
XL	25 - 27 cm	3111792800345
XXL	27 - 30 cm	3111792800352



Thuasne Sport BOA® stabilizing ankle brace

Severe sprain, tendinopathy

035102

BOA®
THE BOA
GUARANTEEUK
CA

Support

INDICATIONS

- Strong ankle support.
- Moderate to severe ankle sprain.
- Joint laxity/chronic instability prevention during sport activities resumption.
- Tendinopathy treatment.
- Conservative treatment of Achilles tendon injuries.

MODE OF ACTION

Ankle stabilization: rigid side stays and elastic straps.

Easy fitting and adjustment to the joint: BOA® Fit System.

COMPOSITION

Polyamide, polyester, neoprene, elastane, polyurethane, stainless steel.



S	19 - 22 cm	3111790351214
M	22 - 26 cm	3111790351221
L	26 - 30 cm	3111790351238

Thuasne Sport elastic wrist support

Wrist sprain, tendinopathy

057503

UK
CA

Support

INDICATIONS

- Wrist sprain.
- Wrist tendinopathies.

MODE OF ACTION

Adjustable fastening system.


COMPOSITION

Polyamide, elastane, synthetic rubber (CR).



Universal size

1	14 - 20 cm	3111795750418
---	------------	---------------



Thuasne Sport thumb strapping

Pain or trauma

033204

UK
CA

Support

INDICATIONS

- Thumb osteoarthritis (including rhizarthrosis).
- Immobilisation in case of pain or trauma of the thumb.

MODE OF ACTION

Thumb stabilization via a plastic wire stay.
Restraint thanks to the anatomical envelope made from elastic fabric.
The splint is held in place by a hook-and-loop strap on the wristband.

COMPOSITION

Polyamide, polyester, synthetic rubber (SBR), polyurethane, elastane.



Universal size

S	14 - 15.5 cm	3111793324314
M	16 - 17.5 cm	3111793324321
L	18 - 20 cm	3111793324338

Thuasne Sport strapping wrist band

Wrist tendinopathies

034002



←
This product contains
latex

UK
CA

Support

INDICATIONS

Wrist tendinopathies.

MODE OF ACTION

Custom closure adjustment with hook and loop fastener.

COMPOSITION

Polyamide, elastane, polyester, latex, polyurethane.



Universal size

1	14 - 20 cm	3111793402395
---	------------	---------------

Thuasne Sport elbow support

Tennis elbow, golfer's elbow, elbow tendinitis

033801



Support

INDICATIONS

- Treatment of lateral or medial epicondylitis (tennis elbow or golfer's elbow).
- Prevention of elbow tendinopathy recurrence.
- Support for physical activity resumption.

MODE OF ACTION

Improved proprioception thanks to the elastic knit.
Compression on the joint.

COMPOSITION

Polyamide, elastane, cotton.



S	22 - 24 cm	3111790338116
M	25 - 27 cm	3111790338123
L	28 - 30 cm	3111790338130
XL	31 - 33 cm	3111790338147

Thuasne Sport anti-epicondylitis elbow brace

Tennis elbow, golfer's elbow, elbow tendinitis

033604



INDICATIONS

- Treatment of lateral or medial epicondylitis (tennis elbow or golfer's elbow).
- Prevention of elbow tendinopathy recurrence.
- Support for physical activity resumption.

MODE OF ACTION

Elbow stabilization.
Improved proprioception.
Even compression on the joint.
Removable tightening strap.

COMPOSITION

Textile components: synthetic rubber (SBR), polyamide, elastane.

Strap: polyester, polyoxymethylene.

Insert: silicone.



S	22 - 24 cm	3111790047285
M	25 - 27 cm	3111790047292
L	28 - 30 cm	3111790047308
XL	31 - 33 cm	3111790047315

Thuasne Sport tennis-elbow armband

Tennis elbow, golfer elbow

035202



UK
CA

INDICATIONS

- Treatment of lateral or medial epicondylitis (tennis/golfer's elbow).
- Joint effusions and swelling due to arthritis and arthrosis.

MODE OF ACTION

Epicondyle support thanks to uniform compression.
Improved proprioception.
Adjustable pressure and precise adjustment by means of the tightening strap.

COMPOSITION

Polyamide, polyurethane, polyester, cotton, elastane, silicone.



1	22 - 30 cm	3111790352112
2	30 - 38 cm	3111790352129



Thuasne Sport Shoulder Pad

Trauma, inflammatory conditions, dislocation **057603 / 04**



UK
CA



Support

INDICATIONS

- Shoulder and/or elbow injuries or inflammatory conditions.
- Shoulder subluxation or dislocation prevention (joint instability).
- Shoulder dislocation.

MODE OF ACTION

Shoulder support thanks to the elastic fabric.
Flat seams provide comfort.
Anatomically-shaped to follow the movement of the shoulder.

COMPOSITION

Polyamide, elastane, synthetic rubber (CR).



S	28 - 31 cm	3111795762015
M	32 - 35 cm	3111795762022
L	36 - 40 cm	3111795762039

Thuasne Sport lumbar support belt

Low back pain, disc/lumbar sprain, back pain **036101**



Support

INDICATIONS

- Non-specific low back pain (acute, subacute and chronic phases).
- Sport activities resumption.
- Sprained disc or lumbar sprain.
- Torn lumbar muscles.

MODE OF ACTION


The device relieves low back pain, supporting the lumbar spine.
2 rigid and conformable stays on the dorsal part.
2 flexible stays on the abdominal part

COMPOSITION

Non-textile components: polyamide (hook and loop),
brushed hardened steel (back stays), polyoxymethylene (ventral stays).
Textile components: polyamide, polyurethane, polyester, elastane.



S	60 - 70 cm	3111790361312
M	71 - 86 cm	3111790361329
L	87 - 102 cm	3111790361336
XL	103 - 121 cm	3111790361343
XXL	121 - 140 cm	3111790361350




SPORTS COMPRESSION

NON MEDICAL DEVICES

**GRADUATED
COMPRESSION
INSPIRED
BY OUR MEDICAL
KNOW-HOW**

Contributing to
correct muscle oxygenation,
slowing down the onset
of cramps and aiding
recovery.





Our products are designed and manufactured in France in our Saint-Étienne factories. Our goal is to improve sports practice and recovery.

Benefits of compression

1.

A MORE EFFICIENT STRIDE
AND REDUCED MUSCLE VIBRATIONS

2.

LESS FATIGUE DURING EXERCISE

3.

PROTECTIVE EFFECT BY REDUCING PAIN

4.

MORE EFFICIENT MUSCLE OXYGENATION

1. Bringard et al., Compression élastique externe et fonction musculaire chez l'homme, Science & Sports (2007), doi: 10.1016/j.scispo.2006.11.003.

2. Clinical report Recovery study July 2015 carried out by THUASNE LPE Univ Savoie Mont Blanc.

3. Clinical report Recovery study July 2015 carried out by THUASNE LPE Univ Savoie Mont Blanc.

4. Clinical report - Recovery July 2015 carried out by THUASNE - LPE - Univ Savoie Mont Blanc.

What is graduated compression?

The most significant pressure is at the ankle and decreases as it rises to below the knee. Thus contributing to good muscle oxygenation, slowing down the appearance of cramps and allowing you to link training sessions more easily.

Good to know



The combination of compression during and after exercise allows you to better recover and link exercise over time.

UP ACTIV Compression high socks

034401



MADE IN FRANCE

ANATOMICAL COMFORT

- Adapted to the morphology of most feet.
- Reinforced mesh on the most stressed areas: heel and toe.
- Reduced friction and slippage in the shoe.
- Limitation of humidity, heat and bad smells.

INNOVATIVE YARNS AND FIBERS FOR OPTIMAL PERFORMANCE

- Helps wick away moisture and reduces foot heating with NILIT® Aquarius weft yarns.
- Reduced bad odors and itching with recycled Polygiène® antibacterial and anti-odor yarns.
- Durable equipment: superior quality fibers that limit wear due to repeated washes.

SUITABLE FOR MOST PEOPLE

UP ACTIV compression high socks are suitable for most levels and morphologies thanks to their targeted mesh and the variety of sizes available.

COMPOSITION


Elastane, polyamide, polyester.


Colors: **Dark blue or sky blue**

Available in 4 sizes: **S, M, L and XL**, depending on the circumference of your ankle

Shoe sizes: **35-38, 39-42 and 43-46**

2 leg heights (normal < 44 cm and long > 44 cm).

Dark blue		35 - 38	39 - 42	43 - 46
S - Normal	19 - 22 cm	3111790344001	3111790344209	-
S - Long	19 - 22 cm	-	3111790344018	3111790344025
M - Normal	22 - 25 cm	3111790344056	3111790344216	-
M - Long	22 - 25 cm	-	3111790344063	3111790344070
L - Normal	25 - 28 cm	3111790344100	3111790344223	-
L - Long	25 - 28 cm	-	3111790344117	3111790344124
XL - Normal	28 - 31 cm	3111790344155	3111790344230	-
XL - Long	28 - 31 cm	-	3111790344162	3111790344179

Sky blue		35 - 38	39 - 42
S - Normal	19 - 22 cm	3111790344032	3111790344049
M - Normal	22 - 25 cm	3111790344087	3111790344094
L - Normal	25 - 28 cm	3111790344131	3111790344148
XL - Normal	28 - 31 cm	3111790344186	3111790344193

UP ACTIV Compression mid socks

034501



MADE IN FRANCE

TECHNICAL CHARACTERISTICS

- Specific mesh tightened at the ankle, allowing degressive compression.
- Reduced friction and slippage in the shoe.
- Helps wick away moisture with NILIT® Aquarius weft yarn.
- Reduced bad odors with recycled anti-bacterial Polygiène® yarn.
- Durable equipment: superior quality fibers that limit wear due to repeated washes.

USE

- Put on the sports compression sock, making sure that the toes and heel are correctly positioned. Unroll the sock to under the calf.


COMPOSITION


Elastane, polyamide, polyester.

Colors: **Dark blue or sky blue**

Available in 4 sizes: **S, M, L and XL**, depending on the circumference of your ankle

Shoe sizes: **35-38, 39-42 and 43-46**

Dark blue		35 - 38	39 - 42	43 - 46
S	19 - 22 cm	3111790345008	3111790345015	3111790345022
M	22 - 25 cm	3111790345053	3111790345060	3111790345077
L	25 - 28 cm	3111790345107	3111790345114	3111790345121
XL	28 - 31 cm	3111790345152	3111790345169	3111790345176

Sky blue		35 - 38	39 - 42
S	19 - 22 cm	3111790345039	3111790345046
M	22 - 25 cm	3111790345084	3111790345091
L	25 - 28 cm	3111790345138	3111790345145
XL	28 - 31 cm	3111790345183	3111790345190

UP ACTIV Compression calf sleeves

034301

RUNTRAIL UPACTIV



MADE IN FRANCE

TECHNICAL CHARACTERISTICS

- Specific mesh tightened at the ankle, allowing degressive compression.
- Reduced friction and slippage in the shoe.
- Helps wick away moisture with NILIT® Aquarius weft yarn.
- Reduced bad odors with recycled anti-bacterial Polygiène® yarn.
- Durable equipment: superior quality fibers that limit wear due to repeated washes.

USE

- The colored border (narrowest part of the sleeve) is placed at the level of the ankle and the wider part under the knee. Put the compression sleeve on the leg. The Thusne logo indicates the front of the sleeve and should be centered below the knee.

COMPOSITION

Elastane, polyamide, polyester.

Colors: **Dark blue or sky blue**Available in 4 sizes: **S, M, L and XL**, depending on the circumference of your ankleShoe sizes: **35-38, 39-42 and 43-46****2 leg heights** (normal < 44 cm and long > 44 cm).

		S	M	L	XL	
Dark blue	Normal	19 - 22 cm	3111790343004	3111790343042	3111790343080	3111790343127
	Long	19 - 22 cm	3111790343011	3111790343059	3111790343097	3111790343134

		S	M	L	XL	
Sky blue	Normal	19 - 22 cm	3111790343028	3111790343066	3111790343103	3111790343141
	Long	19 - 22 cm	3111790343035	3111790343073	3111790343110	3111790343158

UP RECOVERY Long recovery socks

034601

UP RECOVERY



MADE IN FRANCE

TECHNICAL CHARACTERISTICS

- Specific mesh tightened at the ankle, allowing degressive compression.
- Helps wick away sweat with NILIT® Aquarius yarn.
- Reduced odors with Polygiène® recycled anti-bacterial yarn.
- Durable equipment: superior quality fibers that limit wear due to repeated washes.
- Reduced friction and slippage in the shoe.

USE

- Put the compression sock on the foot, making sure that the toes and heel are correctly positioned. Unroll the sock on the leg up to below the knee. Massage the leg to distribute the sock evenly and pull on the tip of the sock to free the toes.

COMPOSITION

Elastane, polyamide, polyester.

Color: **Black**Available in 4 sizes: **S, M, L and XL**, depending on the circumference of your ankleShoe sizes: **35-38, 39-42 and 43-46****2 leg heights** (normal < 44 cm and long > 44 cm).

		35 - 38	39 - 42	43 - 46
S - Normal	19 - 22 cm	3111790346036	3111790346241	-
S - Long	19 - 22 cm	-	3111790346043	3111790346050
M - Normal	22 - 25 cm	3111790346098	3111790346258	-
M - Long	22 - 25 cm	-	3111790346104	3111790346111
L - Normal	25 - 28 cm	3111790346159	3111790346265	-
L - Long	25 - 28 cm	-	3111790346166	3111790346173
XL - Normal	28 - 31 cm	3111790346210	3111790346272	-
XL - Long	28 - 31 cm	-	3111790346227	3111790346234

A HIGH-PERFORMANCE

**ARMSLEEVE,
DYNAMIC COMPRESSION,
INTEGRATED SWEATBAND,
GUARANTEED COMFORT**



Basket ball



Tennis



Handball



Performance armsleeves

034901

Armsleeve

Available
in 2024.



TECHNICAL CHARACTERISTICS

- **The armsleeve combines compression, grip and thickness.**
- **Dynamic compression:** a degressive compression that promotes recovery, proprioception and limits cramps and aches. Less fatigue during exercise.
- **Comfort/Protection:** targeted padded areas: terrycloth wristband with integrated sponge headband. Light vented mesh.
- **Chevron for biceps support.**
- **Limits unpleasant smells:** anti-bacterial fiber.

USE

- **Basketball / Tennis / Handball**

COMPOSITION

Polyamide, elastane.

Colors: **Black and White**

Available in 3 sizes: **S, M, L**

White



	39 - 44 cm	44 - 50 cm	50 - 57 cm
S - Normal : 39-44 cm	14 - 16 cm	3111790348146	-
S - Long : 44-50 cm	14 - 16 cm	-	3111790348153
M - Normal : 39-44 cm	16.5 - 18 cm	3111790348184	-
M - Long : 44-50 cm	16.5 - 18 cm	-	3111790348191
L - Normal : 44-50 cm	18.5 - 20 cm	-	3111790348221
L - Long : 50-57 cm	18.5 - 20 cm	-	-
			3111790348238

Black



	39 - 44 cm	44 - 50 cm	50 - 57 cm
S - Normal : 39-44 cm	14 - 16 cm	3111790348160	-
S - Long : 44-50 cm	14 - 16 cm	-	3111790348177
M - Normal : 39-44 cm	16.5 - 18 cm	3111790348207	-
M - Long : 44-50 cm	16.5 - 18 cm	-	3111790348214
L - Normal : 44-50 cm	18.5 - 20 cm	-	3111790349006
L - Long : 50-57 cm	18.5 - 20 cm	-	-
			3111790349013

Grip Performance long socks

Long socks with non-slip grip

034701

Available
in 2024.



TECHNICAL CHARACTERISTICS

- **The sock combines compression, grip and thickness.**
- **Dynamic compression:** a degressive compression that promotes recovery, proprioception and limits cramps and aches. Less fatigue during exercise.
- **Grip:** non-slip silicone to prevent slippage.
- **Comfort/Protection:** targeted padded areas: terrycloth heel, ankle and toe.
- **Light vented mesh.**
- **Limits unpleasant smells:** anti-bacterial fiber.

USE

- **Basketball / Soccer / Handball / Tennis**

COMPOSITION

Polyamide, elastane, silicone.

Colors: **Black and White**

Available in 6 sizes: **S, S+, M, M+, L, L+**

White				White			
Height < 1.75 m				Height > 1.75 m			
		35 - 38	39 - 42			39 - 42	43 - 46
S	19 - 22 cm	3111790347002	3111790347019	S+	19 - 22 cm	3111790347101	3111790347118
M	22.5 - 25 cm	3111790347040	3111790347057	M+	22.5 - 25 cm	3111790347149	3111790347156
L	25.5 - 28 cm	-	3111790347088	L+	25.5 - 28 cm	3111790347187	3111790347194

Black				Black			
Height < 1.75 m				Height > 1.75 m			
		35 - 38	39 - 42			39 - 42	43 - 46
S	19 - 22 cm	3111790347026	3111790347033	S+	19 - 22 cm	3111790347125	3111790347132
M	22.5 - 25 cm	3111790347064	3111790347071	M+	22.5 - 25 cm	3111790347163	3111790347170
L	25.5 - 28 cm	-	3111790347095	L+	25.5 - 28 cm	3111790347200	3111790347217

Grip Performance mid socks

Mid-calf socks with non-slip grip

034801

Available
in 2024.



TECHNICAL CHARACTERISTICS

- **The sock combines compression, grip and thickness.**
- **Dynamic compression:** a degressive compression that promotes recovery, proprioception and limits cramps and aches. Less fatigue during exercise.
- **Grip:** non-slip silicone to prevent slippage.
- **Comfort/Protection:** targeted padded areas: terrycloth heel, ankle and toe.
- **Light vented mesh.**
- **Limits unpleasant smells:** anti-bacterial fiber.

USE

- **Handball / Soccer / Basketball / Tennis**

COMPOSITION

Polyamide, elastane, silicone.

Colors: **Black and White**

Available in 3 sizes: **S, M, L**

White				White			
		35 - 38	39 - 42			39 - 42	43 - 46
S	19 - 22 cm	3111790347224	3111790347231			3111790348023	3111790348030
M	22.5 - 25 cm	3111790348047	3111790348054			3111790348085	3111790348092
L	25.5 - 28 cm	-	3111790348108			3111790348122	3111790348139

Black				Black			
		35 - 38	39 - 42			39 - 42	43 - 46
S	19 - 22 cm	3111790348016	3111790348023			3111790348023	3111790348030
M	22.5 - 25 cm	3111790348078	3111790348085			3111790348085	3111790348092
L	25.5 - 28 cm	-	3111790348122			3111790348122	3111790348139

CARE & PROTECTION

Thuasne offers a care and protection range to continue to support you in your daily life and your sports activities.



Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

All the medical devices mentioned on this document are CE marked according to the Regulation 2017/745 on medical devices, unless specifically identified as "NOT CE MARKED".

Some of the products featured in this document are also UKCA-marked according to S.I 2002/618 regulation on medical devices.

Devices mentioned in this document are CE class I. Please contact Thuasne should you need any additional information on devices classification.

Please read carefully the instructions for use, indications and contraindications of the product.

**UK
CA**

UK Responsible Person (UKRP):

Thuasne UK Ltd
Unit 4 Orchard Business Centre
North Farm Road
Tunbridge Wells TN2 3XF,
United Kingdom



Thuasne
120 rue Marius AUFAN 92300 Levallois - Perret
www.thuasne.com

Thuasne finger support braces

253501



INDICATIONS

The Thuasne finger support braces provide flexible finger support when returning to sports after injury: severe sprain, fracture.

MODE OF ACTION

- Dynamic buddy splinting: flexible finger support thanks to the 2 independent bands.
- Suitable for long fingers, not suitable for thumb.
- Modular: adaptable tightening and the bands can be cut.
- Interdigital separation to avoid friction between the fingers.
- Non-sticky, flexible and elastic.
- Durable: washable and reusable.

Color: **Navy blue**

COMPOSITION

Polyamide, EVA foam, polyethylene, elastane.
Inside of cup pouch: polyester and Spandex®.
Outside of cup pouch: terylene and cotton.

UNIQUE 3111790056690

Thuasne Sport

Sport Strapping Bandage

048502



INDICATIONS

The Thuasne Sport strapping bandage maintains the muscle or joint weakened by joint, ligament or muscle trauma.

MODE OF ACTION

- Adhesive compression bandage.
- Aerated elastic fabric.
- Multi-directional herringbone coating.
- Short-stretch < 70%.
- Intermediate protective paper.

Color: **White**

COMPOSITION

Cotton, polyester, adhesive mass (synthetic elastomer and tackifying resin), paper.

Width	Length 2.5 m
3 cm	3111790072485
6 cm	3111790072492



**THUASNE SPORT
STRAPPING BANDAGE**

Argicalm®

Reusable hot/cold clay pack

720004

MAIN INDICATIONS**Use cold:**

Sprains, strains, tendinitis, migraines/headache, eye pain, toothache/dental pain, bruises, arthritis, maxillofacial surgery and eye surgery.

Use hot:

Back pain, contractures/muscle and articular stiffness, cramps, torticollis (stiff neck), abdominal pain.

MAIN FEATURES

Practical: reusable hot/cold clay pack. **Easy fitting:** self-fastening elastic strap. **Comfort:** textile protective sleeve. **Anatomical fit:** malleable clay pack. Available in 4 different sizes.

Eye mask and Kids size also available.

SIZE	RECOMMENDED USES	DIMENSION (cm)
Size 1	Wrist - Elbow - Face	19 x 11
Size 2	Ankle - Elbow - Knee	26 x 13
Size 3	Knee - Back - Abdomen	30 x 20
Size 4	Neck - Shoulder - Back	43 x 17
EYES	Eyes	18 x 75
KIDS	Face - Neck - Shoulder - Elbow - Wrist Abdomen - Knee - Ankle	12 x 80

Relieves pain**Hot pain relief****Cold pain relief****GenuStart® Ice**

Cryotherapy knee brace

238901

MAIN INDICATIONS

- Cold and compression therapy following knee ligament injuries and/or ruptures (cruciate and/or lateral ligaments).
- Pre and/or post-operative treatment.
- Knee pain and/or swelling management.

MAIN FEATURES


Cooling effect: 2 sets of 2 reusable ice-packs ensure up to 20 minutes of cooling⁽¹⁾.

Compression: inflatable aircells with external pump.

Ease of application: 2 adjustable self-fastening straps.

Comfort: removable hinged uprights and built-in textile pockets to insert ice-packs.

Bilateral model.



Universal size
35 - 70 cm

(1) Internal CE marking data

Presentation video available online

**Tip!**

The RICE protocol (Rest, Ice, Compression and Elevation) is recommended for treating inflammation after trauma or surgery. It is the preferred treatment for the first 4 to 5 days post-op.⁽¹⁾

...nce for rest, ice,
...nt of ankle sprains in

**GenuStart® Ice**

Thuasne Sport Reseller Opportunities

Become a Thuasne Sport reseller to access competitive discounts with no minimum order quantity, credit terms, exclusive access to branding and marketing assets and an easy to use online ordering portal.

Contact our Customer Service Team to enquire about your seller account today.

e: orders@thuasne.co.uk
t: 01892 481620

SPORT

Elastic Wrist Support

Reinforced Knitted Knee Support

WINGS FOR YOUR HEALTH

The Thuasne Sport range can be viewed and bought online from:
www.thuasne-store.co.uk



For technical advice on our products or to place an order, contact
Thuasne Customer Service:
Telephone: 01892 481620
Email: orders@thuasne.co.uk



Head office: Thuasne UK Ltd., Unit 4 Orchard Business Centre, North
Farm Road, Tunbridge Wells, Kent, TN2 3XF

www.thuasne.co.uk